



VIVID DIGITAL DENTURES

3X Stronger than traditional
Multi-Layer milled teeth
Superior quality

DIGITAL RECORD ON FILE

Copies can be made up to 3
years from the date of the
original case without having to
bring the patient back in

VIVID DIGITAL DENTURE ECONOMY

- Economy Multi Layered milled teeth
- Lucitone Print Denture Base

VIVID DIGITAL DENTURE PREMIUM

- Premium Multi Layered milled teeth
- Lucitone Print Denture Base

VIVID DIGITAL DENTURE ULTRA

- Premium Multi Layered milled teeth
- Ivoclar Milled Denture Base

Replacing a current denture or prescribing a new denture:

Step One: Impression/Bite

Step Two: Tryin (optional)- Monolithic printed

Step Three: Final

What to send: Scan of existing denture with wash
reline, bite scan, opposing scan, and photos of
smile if possible.

Immediate denture:

Step One: Scan of arch with teeth

Step Two: Tryin (optional)- Monolithic printed

Step Three: Final

What to send: Scan of existing dentition, bite scan,
opposing scan, and photos of smile if possible.

Copy of current denture:

Step One: Scan of current denture/ or send
patient to the lab with their denture

Step Two: New Copy Denture delivered

What to send: Scan of existing denture. With wash
reline, if ill fitting