



GPS
DIGITAL DENTAL LAB

Since 2019, GPS Digital Dental Lab has embraced digital denture technology, using Digital Print Resins and Milling Methods for all IOS-scanned cases. The results have been outstanding—patients enjoy better fit, enhanced esthetics, and superior strength compared to traditional acrylic dentures. Plus, digital design ensures easy reproduction if lost or damaged.

Given this success, we are transitioning our entire removable department to fully digital dentures. We're confident that you and your patients will love the results.

If you use an IOS scanner, you may already be familiar with our digital dentures. For best practices, scan the QR codes below for tutorial videos.

Reference Denture



Copy Denture



Immediate Denture



For those still using traditional analog methods, no changes are necessary on your part. Once we receive the models articulated from your custom tray and bite rim, we will scan those models and digitize everything to design your full denture, from the **prototype try-in** to the final denture. To see what to expect when you receive your try-in and to learn more about the differences between our VIVID DIGITAL DENTURES, please scan the QR code below. If you would like to receive samples, please reach out to the lab, and we will be happy to accommodate you.

Thank you,

Andrew Perricone

President/CEO

aperricone@gpsdental.com



GPS
DIGITAL DENTAL LAB



3501 Parkway Center Court Orlando, FL 32808

Office: 407-447-4250

www.gpsdental.com